

pure
enrichment

5 YEAR
WARRANTY

pure relief™

CORDLESS LUMBAR & ABDOMINAL
HEATING WRAP



USER MANUAL
MODEL: PECLPBK-MG

IMPORTANT SAFETY INSTRUCTIONS DO NOT DESTROY

DANGER: TO REDUCE THE RISK OF BURNS, ELECTRIC SHOCK, FIRE, AND PERSONAL INJURY, THIS PRODUCT MUST BE USED IN ACCORDANCE WITH THE FOLLOWING INSTRUCTIONS:

1. READ ALL INSTRUCTIONS.
2. DO NOT USE WHILE SLEEPING.
3. DO NOT USE ON AN INFANT OR ON AN ANIMAL.
4. THIS PRODUCT IS NOT TO BE USED ON OR BY AN INVALID, A SLEEPING OR UNCONSCIOUS PERSON, A PERSON WITH DIABETES, OR A PERSON WITH POOR BLOOD CIRCULATION UNLESS CAREFULLY ATTENDED.
5. DO NOT USE ON AREAS OF THE SKIN THAT ARE INSENSITIVE TO TOUCH OR HEAT.
6. BURNS CAN OCCUR REGARDLESS OF CONTROL SETTING. CHECK SKIN UNDER PRODUCT FREQUENTLY.
7. DO NOT USE IN AN OXYGEN-FILLED ATMOSPHERE.
8. DO NOT USE PINS OR OTHER METALLIC OBJECTS TO FASTEN THIS PRODUCT INTO PLACE.
9. DO NOT SIT ON OR AGAINST THE PRODUCT OR CRUSH IT. AVOID SHARP FOLDS IN THE FABRIC.
10. DO NOT USE THE CORD AS A HANDLE.
11. CHILDREN SHOULD NOT USE THIS UNIT WITHOUT ADULT SUPERVISION.

12. PLACE UNIT ON TOP OF AND NOT UNDER THE PART OF THE BODY THAT REQUIRES HEAT THERAPY.
 13. ONLY USE THIS PRODUCT ON A 110-120V AC CIRCUIT.
 14. UNPLUG WHEN NOT IN USE.
 15. DO NOT USE WITH GENERATORS, POWER CONVERTERS, OR INVERTERS.
 16. DO NOT TAMPER WITH THIS PRODUCT IN ANY WAY.
 17. THERE ARE NO USER-SERVICEABLE PARTS. IF FOR ANY REASON THIS PRODUCT DOES NOT FUNCTION PROPERLY, CONTACT PURE ENRICHMENT AT HELP@PUREENRICHMENT.COM
 18. DO NOT USE THIS PRODUCT WITH LINIMENT, SALVE OR OINTMENT PREPARATIONS AS SKIN BURNS CAN RESULT.
 19. NEVER PULL THIS PAD BY THE SUPPLY CORD.
 20. DO NOT BEND OR PINCH THE CORD.
 21. DO NOT USE PRODUCT IF THE CORD IS DAMAGED.
 22. NEVER LEAVE THE PRODUCT UNATTENDED, ESPECIALLY IF CHILDREN ARE PRESENT.
 23. LOOP CORD LOOSELY WHEN STORING. TIGHT WRAPPING CAN DAMAGE CORD & INTERNAL PARTS.
 24. THIS PRODUCT HAS A POLARIZED PLUG (ONE BLADE IS WIDER THAN THE OTHER). TO REDUCE THE RISK OF ELECTRIC SHOCK, THIS PLUG IS INTENDED TO FIT IN A POLARIZED OUTLET ONLY ONE WAY. IF THE PLUG DOES NOT FIT FULLY IN THE OUTLET, REVERSE THE PLUG. IF IT STILL DOES NOT FIT, CONTACT A QUALIFIED ELECTRICIAN TO INSTALL THE PROPER OUTLET.
- DO NOT ATTEMPT TO MODIFY THE PLUG OR DISABLE THIS SAFETY FEATURE IN ANY WAY.
25. SAVE THESE INSTRUCTIONS.

FCC NOTICE

Model: PureRelief™ Cordless Lumbar & Abdominal Heating Pad (PECLPBK-MG)

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Pure Enrichment
2803 S Yale St.
Santa Ana, CA 92704
pureenrichment.com

The product is protected against unauthorized use. DO NOT modify it without the authorization of the manufacturer, or the user's authority to operate could be voided.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase separation between equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

INTRODUCTION

Thank you for purchasing the PureRelief™ Cordless Lumbar & Abdominal Heating Wrap from Pure Enrichment. This specialized wrap can be used on the go to provide temporary targeted heat therapy to the lower back, abdomen, knee, arm, and shoulders.

PRODUCT FEATURES

We designed the PureRelief™ Cordless Lumbar & Abdominal Heating Wrap with several features that make it the ideal solution for comfortable daily use. These features include:



CORDLESS USE

Rechargeable battery pack allows you to use the heating pad at school, work, and on-the-go!



4 HEAT SETTINGS

Choose from warm, low, medium, or high heat.



OPTIONAL GEL PACK

Includes a therapeutic gel pack that can be warmed or cooled and inserted into the built-in gel pack pocket.



SUPER SOFT

Luxurious and plush micromink fabric.



2-HOUR AUTO SHUT-OFF

Shuts off heat automatically after 2 hours of continuous use.

CHARGING PureRelief™

1. Attach the DC plug of the adapter to the Battery Pack. Push both parts together securely to ensure a tight fit.
2. Plug the adapter into a 120V AC electrical outlet.
3. While charging, the 4th heat level LED light will flash on and off in red. When fully charged (after about 2.5 hours) the 1st heat level LED light will turn on and stay green.

CHARGING PureRelief™ (CONT'D)

- After charging, disconnect adapter from the Battery Pack.

NOTE: If the heat level indicators do not light up or stay on, make sure the DC plug of the adapter is securely plugged into the Battery Pack.

USING PureRelief™



NOTE: Fully charge the Battery Pack before first use.

The Battery Pack must be connected to the Heating Wrap Connector before the power will turn on.

- Attach the Battery Pack to the Heating Wrap Connector in the Battery Pack Pocket. Push both parts together securely to ensure a tight fit. Place the Battery Pack into the Pocket.
- To turn the unit on, press the power button until the 1st heat level LED indicator light turns on.
- Select the desired heat setting (1-warm, 2-low, 3-medium, or 4-high) by continuing to press the power button until the heat indicator light turns on your desired setting.
- To turn off, press the power button until the lights turn off. Wrap turns off automatically after 2 hours of continuous use.

NOTE: If the heat level indicators do not light up or stay on, make sure the Battery Pack is securely plugged into the Heating Pad Connector.

INSTALLING HOT/COLD GEL PACK

PureRelief™ is made with a built-in pocket that holds a hot/cold gel pack. For cold use, freeze the gel pack for 2 hours before the first use and for 20 minutes thereafter. For warm use, microwave the gel pack for 1 minute. If it is not warm enough, microwave in 10-second intervals until the desired temperature is achieved.

APPLYING PureRelief™

PureRelief™ can be used on your lower back, abdomen, arm, shoulder, or knee. To apply, wrap the pad around the desired treatment area making sure the gel pack side is touching the treatment area. Secure the pad in place by firmly pressing on the fasteners. Do not fasten too tightly.

WASHING INSTRUCTIONS

DO NOT DRY CLEAN! DO NOT USE BLEACH! DO NOT USE WRINGER!

- Disconnect PureRelief™ from the Battery Pack.
- Machine wash PureRelief™ in cold water on the gentle cycle, or hand wash. Tumble dry on low heat or line dry.
- Do not plug in PureRelief™ until the Heating Pad Connector is completely dry.

PRODUCT SPECIFICATIONS

Size:	51.2" x 9.5"
Weight:	Heating Wrap: 4 lb, Gel Pack: 13 oz
Power Supply:	AC: 100-240V, 50/60Hz, 15W±10%
Pad Material:	Micromink [100% polyester (both sides)]
4 Heat Settings:	Warm—green light, temperature range (~110 °F) Low—yellow light, temperature range (~120 °F) Medium—orange light, temperature range (~130 °F) High—red light, temperature range (~140 °F)
Battery Usage Time:	1—Warm—green light ~140 minutes 2—Low—yellow light ~110 minutes 3—Medium—orange light ~80 minutes 4—High—red light ~60 minutes

WARRANTY

PureRelief™ comes with an industry-leading 5 Year Warranty that begins on the date of purchase.

The warranty applies to the heating pad and necessary parts and labor relating thereto. The warranty does not apply to damage resulting from failure to follow the operating instructions, accidents, abuse, alterations, or disassembly by unauthorized individuals.

PRODUCT REGISTRATION

To ensure your full warranty and to receive product updates and streamlined customer support, remember to register your product at:

pureenrichment.com/productregistration

ARE YOU 100% SATISFIED?

If you have product questions or concerns, don't hesitate to contact us directly at:

Email: help@pureenrichment.com

Phone: (657) 275-3737

(Available Monday-Friday, 8am-5pm PST)

Pure Enrichment has an award-winning warranty, exchange, and customer service program that guarantees hassle-free solutions to any issue you might have within 24 hours!



Santa Ana, CA 92704
pureenrichment.com