



USER'S GUIDE

Smart Bluetooth Baby Scale

Hubble Grow

The Features described in this user's guide are subject to modifications without prior notice.

Welcome...

to your new Hubble Smart Baby Scale!

Thank you for purchasing your new Hubble Smart Baby Scale.

The digital baby scale has been conceived to control the baby's growth and quantify the milk taken at each feed. With a maximum load of 30 kg and a division of 10 g, the scale allows to follow the evolution of the baby's weight during his growth.

The pack comes with a soft pad for placing on the scale so that it gives comfort to your love one when weighing.

You can also use our Hubble Connected app to track your baby's growth, it will also calculate baby's BMI according to input weight & height.

Please keep your original dated sales receipt for your records. For warranty service of your Hubble product, you will need to provide a copy of your dated sales receipt to confirm warranty status. Registration is not required for warranty coverage.

For product related questions, please contact:

Live Chat: hubbleconnected.com

Email: support@hubbleconnected.com

<https://hubbleconnected.com/warranty>


<https://hubbleconnected.com/apps/product-registration>

Visit **support.hubbleconnected.com** to download the User Guide with instructions, answers to frequently asked questions, and more to help you get the most from your product.

Inside the Box

- 1 x Baby Scale
- 1 x Quick Start Guide
- 1 x Soft Pad

Overview of the control panel



1. LCD Display

Bluetooth Signal

Memory symbol

Flashes when the memory is full.

Weight Change indicator

The symbol appears together with the Memory symbol to show weight value change with Up arrow (weight gain) or Down arrow (weight loss) indications.

2. Power ON/OFF button

3. UNIT button

Select the unit of measurement: kilograms (kg) or pounds (lb).

4. TARE button

Reset the scale to zero reading.

5. HOLD button

If the baby is moving, press the Hold button to lock in the correct weight. Once you pressed, the data will be stored.

Press HOLD while the baby is on the platform alone and still, so the machine stabilises the detected weight (it will flash on the screen). When the display stops flashing, the weight of the baby will be displayed alternating with the up/down difference from the last weight.

Table of contents

1. Safety Instructions	5
2. Getting Started	7
2.1 INSTALL BATTERIES	7
3. Using the Baby Scale	8
3.1 Turn On/Off the baby scale	8
3.2 Select Weighing Unit	8
3.3 Tare Weight	8
3.4 Weight Tracker	8
3.4.1 First weighing	8
3.4.2 Following weighings	8
4. Troubleshooting	10
5. Disposal of the Device (Environment)	11
6. Cleaning	12
7. General Information	13
8. Technical Specifications	17



1. Safety Instructions

WARNING:

This product is compliant with all relevant electromagnetic field standards and is, when handled as described in the User's Guide, safe to use. Therefore, always read the instructions in this User's Guide carefully before using the device.

- Do not drop the Scale.
- Do not jump on the scale.
- Do not weigh loads totaling more than 66lbs (30kgs).
- Do not expose the Scale to excess moisture or humidity.
- Place the weighing scale on a completely flat and rigid surface. A soft and uneven surface (e.g. carpets or irregular surfaces) can compromise the scale's precision.
- Careful that all 4 feet rest on the flat surface since they are equipped with load sensors.
- This product is not a toy. Do not allow children to play with it.
- Keep this User's Guide for future reference.
- Test the product and all its functions so that you are familiar with it prior to use.
- Do not install the Scale near a heat source.
- Clean the platform with a wet cloth, clean the bottom of the scale with a dry and soft cloth.
- Do not use chemical or abrasive products.
- The machine is not waterproof. Do not let liquids enter in the scale's body. Do not leave or use the scale in waterlogged environments. Humidity over 85% or contact with water or other liquids can cause corrosion and compromise the scale's correct operating.



- Only use the same battery type This product uses replaceable battery AAA sized battery cells (LR03).
- Do not touch the battery contacts with sharp or metal objects.
- Cleaning and user maintenance shall not be made by children without supervision. Cleaning and user maintenance shall not be made by children without supervision.

CAUTION


There is danger of explosion if the battery is replaced by an incorrect type. Dispose of used batteries according to their instructions.



2. Getting Started

2.1 Install Batteries

The battery operates with 4 AAA alkaline batteries type LR03 (1.5V), not included.

1. The battery housing is on the bottom of the weighing scale, carefully pry open the battery cover on the underside of the Scale.
2. Install four (4) new "AAA" batteries (LR03) into the battery compartment, making sure to match the positive and negative ends of the batteries to the correct contacts as indicated in the battery compartment.
3. Replace the battery cover after insert the batteries.
4. Press the Power ON/OFF  button to turn on the weighing scale.


Note:

- Batteries are to be installed with the correct polarity.
- Do not mix batteries of different types or mix old and new batteries, replace all 4 batteries when empty with new batteries.
- If the scale is to be stored unused for a long period of time, the batteries should be removed.
- The battery contacts are not to be short-circuited.




3. Using the Baby Scale

3.1 Turn On/Off the baby scale


Press the Power ON/OFF  button to turn On/Off the weighing scale.

3.2 Select Weighing Unit

Press the UNIT  button on the control panel to select unit of measurement, kilograms (kg) or pounds (lb).

3.3 Tare Weight



The weighing scale zeros automatically when turned on and displays the digit "0.00", if you weight the baby with a soft cloth or towel, it need to tare weight.

1. To add a towel or cloth when the machine is already turned on.
2. Press the TARE  button to zero the weight of soft cloth or towel placed on the platform, and weight that the digit "0.00" is displayed.


3.4 Weight Tracker

The digital baby scale allows to follow the evolution of the baby's weight during his growth, with a maximum load of 30 kg and a division of 10 g.

3.4.1 First weighing

1. Press Power ON/OFF  button to turn on the baby scale, when the digits "0.00" are displayed, the machine is ready for weighing (If necessary, use the function TARE as previously described section 4.3).
2. Place the baby on the weighing scale and the weight will be displayed on the screen.
3. Press the HOLD  button to stabilize the readings , and the weight data will be saved in the memory of baby scale.

3.4.2 Following weighings

1. To turn on the baby scale by pressing Power ON/OFF  button.
2. Place the baby on the weighing scale when the digit "0.00" is displayed.




3. The weight of the baby will appear on the screen, press the HOLD  button to stabilize the readings.
4. The measured weights stop flashing and following weight data will display as alternating.
 - Weight of baby
 - Weight Change (Compare with last weighing)



Fig.1

- if it is weight loss, the symbol ↓ and weight change are display on the screen (Fig.1).



Fig.2

- if it is weight gain, the symbol ↑ and weight change are display on the screen (Fig.2).



Fig.3

- No weight change was detected (Fig.3).

6. The weight data will be saved in the memory of baby scale, the weighing scale with turn off after 20 seconds if no any operation.

NOTE:

The weight data saved in scale memory will be erased when batteries are removed or replaced, the weighing scale will start again from "First Weighing" to record the weight data.



4. Troubleshooting

When a malfunction or incorrect weighting measurement occurs, an error message will appear as described below.

Error Code	Problem Description	Solution
Lo	Low battery.	Replace all 4 batteries with new ones.
Err	The weight exceeds the scale's maximum load.	Please make sure not to exceed the maximum load of 30 kg.
Err	Wrong operating.	Read through the User's Guide or Quick Start Guide, and follow the instructions and Warnings when operating.
No Display	No Power Display time out	<ul style="list-style-type: none"> • Try replacing the AAA battery cells. • Check if the Scale is turned on. Add weight on the Scale to turn it ON. • The display will go blank in 8s after you remove the weight from the scale.
Inaccurate Measurement	Improper placement	<ul style="list-style-type: none"> • Make sure the Scale is placed on a flat and hard surface, such as hardwood or ceramic tile flooring.



5. Disposal of the Device (Environment)

At the end of the product life cycle, you should not dispose of this product with normal household waste. Take this product to a collection point for the recycling of electrical and electronic equipment. The symbol on the product, user's guide and/or box indicates this obligation.



By taking end-of-life products to a recycling point, you are ensuring that some of the product materials can be re-used. By reusing some parts or raw materials from used products, you make an important contribution to the protection of the environment.

Exhausted batteries do not belong with household waste. Dispose of batteries in accordance with the current federal, state and local regulations. As a consumer, you are legally obligated to return used batteries.

Note:

1. *In environments where there is a high frequency of electrostatic discharge, the unit may malfunction and the user may have to reset it.*
2. *Dispose of used batteries with care; please consult retailer for details. Please contact your local authorities should you require additional information on the collection points in your area.*



6. Cleaning

Clean your Scale with a slightly damp cloth or with an antistatic cloth. Never use cleaning agents or abrasive solvents.

Cleaning and care

- Do not clean any part of the Scale with thinners or other solvents and chemicals - this may cause permanent damage to the product which is not covered by the warranty.
- Keep your Scale away from hot, humid areas or strong sunlight, and do not let it get wet.
- Every effort has been made to ensure high standards of reliability of your Scale. However, if something does go wrong, please do not try to repair it yourself - contact Customer Service for assistance.



7. General Information

If your product is not working properly, read this User Guide.

Contact Customer Service:

Live Chat: hubbleconnected.com

Email: support@hubbleconnected.com

Online warranty registration: <https://hubbleconnected.com/warranty>

Online Product Registration:

<https://hubbleconnected.com/apps/product-registration>

Consumer Products and Accessories Limited Warranty ("Warranty")

Thank you for purchasing this Hubble Smart Nursery branded product manufactured under license by Binatone Electronics International LTD ("BINATONE").

What Does this Limited Warranty Cover?

Subject to the items listed below in the section "What is not Covered" BINATONE warrants that this Hubble connected branded product ("Product") will be free from defects in material and workmanship under normal use for a period of 1 year from the date of purchase and that any certified accessory ("Accessory") sold for use with this product will be free from defects in materials and workmanship under normal use for a period of 90-days from the date of purchase. This Limited Warranty is your exclusive warranty and is not transferable.

What is not covered?

This Limited Warranty excludes and does not cover:

- Normal deterioration or wear resulting from use;
- Damage caused by or resulting from misuse, abuse, or accident;
- Damage caused by or resulting from improper operation or storage;
- Damage caused by or resulting from contact with foreign substances, including, but not limited, to liquid, water, rain, humidity, perspiration, sand, dirt, and/or food;



- Damage caused by or resulting from the use of Non-Hubble Smart Nursery branded or certified Products or Accessories or other peripheral equipment;
- Damage caused by or resulting from modifications, alterations, tampering, improper installation, or repairs by someone other than Hubble connected, BINATONE or their authorized representatives;
- Defects, damages, or the failure of Products or Accessories due to any communication service or signal you may subscribe to or use with the Products or Accessories; and,
- Periodic maintenance, repair and replacement of parts due to normal wear and tear.

What Other Limitations Are There?

THIS LIMITED WARRANTY SETS FORTH THE ENTIRE LIABILITY OF BINATONE AND CONTAINS YOUR SOLE AND EXCLUSIVE REMEDY. OTHER THAN THIS LIMITED WARRANTY, BINATONE MAKES NO ADDITIONAL REPRESENTATIONS OR WARRANTIES, EXPRESS OR IMPLIED, INCLUDING WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IF YOUR STATE DOES NOT PERMIT THE DISCLAIMER OF WARRANTIES, THEN THE IMPLIED WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE IS LIMITED TO THE DURATION OF THIS LIMITED WARRANTY. REPAIR OR REPLACEMENT IS YOUR SOLE AND EXCLUSIVE REMEDY AGAINST BINATONE AND Hubble Smart Nursery. IN NO EVENT SHALL BINATONE OR Hubble Smart Nursery's LIABILITY TO REPAIR OR REPLACE YOUR PRODUCT OR ACCESSORY EXCEED THE ORIGINAL PURCHASE PRICE OF THE PRODUCT AND ACCESSORY AND BINATONE AND Hubble Smart Nursery's TOTAL LIABILITY FOR DAMAGES ARISING FROM THIS LIMITED WARRANTY SHALL NOT EXCEED THE PRICE YOU PAID FOR THE PRODUCT AND ACCESSORY REGARDLESS OF WHETHER SUCH CLAIM ARISES UNDER CONTRACT, STATUTE, TORT OR OTHERWISE. IN NO EVENT AND UNDER NO CIRCUMSTANCES WILL BINATONE OR Hubble Smart Nursery BE LIABLE FOR DIRECT, INDIRECT, INCIDENTAL,



SPECIAL OR CONSEQUENTIAL DAMAGES OF ANY KIND, OR LOSS OF REVENUE OR PROFITS, LOSS OF BUSINESS, LOSS OF INFORMATION OR OTHER FINANCIAL LOSS ARISING OUT OF OR IN CONNECTION WITH THE ABILITY OR INABILITY TO USE THE PRODUCTS OR ACCESSORIES. ALL DISCLAIMERS, LIMITATIONS OF LIABILITY AND MODIFICATIONS OF WARRANTIES STILL APPLY EVEN IF THE LIMITED REMEDY OF REPAIR AND REPLACEMENT FAILS OF ITS ESSENTIAL PURPOSE.

Some jurisdictions do not allow the limitation or exclusion of incidental or consequential damages, so the above limitations or exclusions may not apply to you. This Limited Warranty gives you specific legal rights, and you may have other rights that vary from one jurisdiction to another.

Who is covered?

This Warranty extends only to the first consumer purchaser and is not transferable.

What will BINATONE do?

BINATONE or its authorized distributor at its sole discretion will, at no charge, repair or replace any Products or Accessories that do not conform to this Warranty. In the event that replacement is necessary to remedy a defect in material or workmanship, BINATONE may use new or remanufactured products, accessories or parts of equal quality to the original Product or Accessory.

How to Obtain Warranty Service or Other Information?

To obtain service or information, please contact us by:

Live Chat: hubbleconnected.com

Email: support@hubbleconnected.com

Online warranty registration: <https://hubbleconnected.com/warranty>

Online Product Registration:

<https://hubbleconnected.com/apps/product-registration>

You will receive instructions on how to ship the Products or Accessories at your expense and risk, to a BINATONE Authorized Repair Center.

You must include: (a) the Product or Accessory; (b) the original proof of purchase (receipt) which includes the date, place and seller of the Product;



(c) if a warranty card was included in your box, a completed warranty card showing the serial number of the Product; (d) a written description of the problem; and, (e) your name, address, email address and telephone number. These terms and conditions constitute the complete warranty agreement between you and BINATONE regarding the Products or Accessories purchased by you, and supersede any prior agreement or representations, including representations made in any literature publications or promotional materials issued by BINATONE or representations made by any agent, employee or staff of BINATONE, that may have been made in connection with your purchase of any Product and/or Accessory.

FCC and IC regulations (FCC Part 15)

Note: *This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the Federal Communications Commission (FCC) rules.*

Industry Canada. This device complies with Industry Canada license-exempt RSS standard(s): Operation is subject to the following Two conditions:(1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

Download User Guide from website: support.hubbleconnected.com



8. Technical Specifications

Weigh range	0.22lbs(100g) - 66lbs(30kg)
Measurement Accuracy	0.02lbs (10g)
Power supply	4 x 1.5V AAA size alkaline batteries (LR03) NOT INCLUDED
Operating environment	10°C- 35°C / 50°F- to 95°F with relative humidity up to 85% (non condensing)
Storage/ transportation environment	10°C- 50°C / 50°F- to 122°F with relative humidity up to 95% (non condensing)
Dimensions	approx. 21.22"x13"x3.35" (539x330x85mm) LxWxH
Weight	approx. 3.52lbs (1.6kg)





Manufactured, distributed and sold by Hubble Baby(HK) Limited. Hubble Connected and the logo are trademarks or registered trademarks of Binatone Electronics International LTD. All rights reserved.

US_EN Version 1.0
Printed in China

