



Questions?
800 345 4449
info@cotrancorp.com
cotrancorp.com


COTRAN
COTRAN CORP
45 High Point Ave
Portsmouth, RI
© 2019

You can easily remove the insert for cleaning. **Warning: Do not try to clean the membrane inside the insert. Just dip it in cleaning solution**

Cleaning

The implementation of an insert inside the BESS nipple allowed us to introduce a round milk opening at the tip. This 2mm hole controls the milk intake and substantially increases the life of the nipple by eliminating tearing of the nipple by eliminating tearing of the traditional cross design.

Control Flow



item# 15-310



Screw-on nipple
with insert and ring cap



item# 15-310

Leak-Free • Tear-Free • Long-Lasting