

# K&H Original Bolster Pet Cot™

## Care & Instructions

Made in China

K&H Manufacturing, LLC  
Colorado Springs, CO  
www.khmf.com | (719) 591-6950



**Thank you for purchasing K&H Manufacturing's Original Bolster Pet Cot™**

**1**

Fold out cover so it is flat.

**2**

Be sure all 4 frame tubes are threaded through the sleeves sewn along edges of the cover.

**3**

Insert all 4 frame tubes into the nylon corner supports.

**4**

Line up the bolster on the top of the pet cot. Slide the back legs of the pet cot through the elastic straps. Fasten each hook and loop fastener on the front of the bolster around the metal frame until all are secure and tight.

With proper care, your K&H Original Bolster Pet Cot™ should provide years of comfort and service.

### Care Instructions:

Wipe clean with a damp cloth.  
Machine wash gentle, line dry.



### K&H Limited One Year Warranty

**Our warranty is as follows:** K&H Manufacturing warrants to the original consumer, that this product will be free of defects and workmanship for a period of one year from the date of purchase. The manufacturer's liability hereunder is limited to replacement of the product. This warranty is void if the product has been damaged by accident, unreasonable use, neglect, tampering, or other causes not arising from defects in material or workmanship.

**To exercise your warranty:** If, during the warranty period, your product malfunctions, return it to the place of purchase with your **original purchase receipt**. If this cannot be done, return the product, postage prepaid, to K&H Manufacturing, with **dated proof of purchase or original date receipt**. A return authorization number, obtained by calling (719) 591-6950 must accompany all returns to manufacturer. Please include a note explaining the problem when returning the product.